

TO OPEN

- Remove Lid from carrier.
- Lift flap in center of lid.
- Pull up corner of center wipe, twist to a point and thread through dispenser hole in lid. DO NOT PUSH FINGER THROUGH OPENING.
- Replace lid. Pull out exposed wipe and snap off. The next wipe pops up automatically. Dispose of used wipe in trash receptacle. Snap center flap down when finished to retain moisture.
- Do not flush.

FEATURES

- Ready to use wipes - Just pull, wipe & toss
- Antibacterial formula cleans, disinfects & sanitizes in one step (see chart for all contact times)
- Hospital & General disinfectant wipes
- Made with scrubbing fibers
- Cleans without scratching your surfaces
- Cuts through tough grime
- Use every day in every room throughout the house
- Kills Antibiotic Resistant Bacteria: CA-MRSA, MRSA, VISA, VRE & NDM-1
- Kills Influenza A2/Japan in 2 minutes

DIRECTION FOR USE

CLEANING: Use one *Disinfectant Wipes* to completely remove surface of all major dirt and debris from surfaces to be treated.

DISINFECTING: To disinfect hard, non-porous surfaces, thoroughly wet the surface to be treated. Treated surface must remain visibly wet for disinfecting. Use enough wipes for treated surface to remain visibly wet for 10 minutes or refer to the contact times indicated in the chart. For heavily soiled surfaces, clean before following disinfecting instructions. Avoid contact with food.

Remove all animals/poultry and their feed from premises prior to disinfection. Remove all heavy soil, such as urine and fecal matter, from surfaces and objects prior to disinfection. Empty all feeding and watering appliances prior to disinfection. Following disinfection, do not house animals/poultry until areas have been ventilated. All surfaces and/or objects that will contact feed or drinking water should be rinsed with potable water before reuse. **SANITIZING/DEODORIZING:** To sanitize/deodorize hard, non-porous surfaces; wipe surface, using enough wipes to allow treated surface to remain wet for 1 minute for *Staphylococcus aureus* (Staph) (ATCC 6538), *Klebsiella pneumoniae* (ATCC 43212), *Listeria monocytogenes* (ATCC 35291), *Escherichia coli* (E. coli) O157:H7 (ATCC 35150), *Campylobacter jejuni* (ATCC 29428). Let surface dry. For highly soiled surfaces, clean excess dirt first.

FUNGICIDAL ACTIVITY: *Disinfectant Wipes* is an effective fungicide against *Trichophyton interdigitale* (ATCC 9633), when used as directed on hard surfaces found in bathrooms, shower stalls, locker rooms, exercise facilities or other clean, hard, non-porous, surfaces commonly contacted by bare feet. The surface must remain visibly wet for 3 minutes. Allow surfaces to air dry.

Disinfectant Wipes is a bactericidal in the presence of organic soil (5% blood serum) against *Escherichia coli* (E. coli) O157:H7 (ATCC 35150), *Salmonella enterica* (ATCC 10708), *Pseudomonas aeruginosa* (ATCC 15442), *Staphylococcus aureus* (Staph) (ATCC 6538), *Methicillin resistant Staphylococcus aureus* (MRSA), *Community Associated Methicillin Resistant Staphylococcus aureus* (CA-MRSA) (Genotype USA400), *Community Associated Methicillin Resistant Staphylococcus aureus* CA-MRSA (Genotype USA300), *Klebsiella pneumoniae* New Delhi Metallo-Beta Lactamase (NDM-1) Carbapenem Resistant, *Vancomycin intermediate resistant Staphylococcus aureus* (VISA), at a 2 minute contact time. *Streptococcus pyogenes* (Strep) (ATCC 19615), *Vancomycin resistant Enterococcus faecium* (VRE), at 5 minutes contact time.

Disinfectant Wipes is an effective disinfectant against *Escherichia coli* O157:H7 (ATCC 35150), *Pseudomonas aeruginosa* (ATCC 15442), *Salmonella enterica* (Salmonella) (ATCC 10708), *Staphylococcus aureus* (ATCC 6538) & *Streptococcus pyogenes* (Strep) (ATCC 19615), *Community Associated Methicillin Resistant Staphylococcus Aureus* (CA-MRSA) (Genotype USA400) & *Community Associated Methicillin Resistant Staphylococcus Aureus* (CA-MRSA) (Genotype USA300), *Vancomycin Resistant Enterococcus faecium* (VRE), *Methicillin Resistant Staphylococcus Aureus* (MRSA), *Vancomycin Intermediate resistant Staphylococcus Aureus* (VISA), *Klebsiella pneumoniae* New Delhi Metallo-Beta Lactamase (NDM-1) Carbapenem Resistant.

Disinfectant Wipes is an effective disinfectant on hard, non-porous surfaces against influenza A2/Japan, Herpes Simplex Virus 1 (ATCC VR-260) and Herpes Simplex Virus 2 (ATCC VR-734) when the treated surface is allowed to remain wet for 2 minutes (30 seconds for SARS-CoV-2 (the virus that causes COVID-19)).

KILLS HIV-1, HBV & HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings (Hospitals, Nursing Homes) or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus, Type 1, Hepatitis B Virus (HBV) & Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR DECONTAMINATION & CLEANING AGAINST HIV-1, HBV & HCV OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

PERSONAL PROTECTION: When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks, and eye coverings.

CLEANING PROCEDURES: Target surfaces must be cleaned prior to disinfection, as the presence of heavy soil on a surface may reduce the intended disinfectant efficacy of the product.

CONTACT TIME: Allow surfaces to remain wet for 2 minutes. (HIV-1 is inactivated in 30 seconds, HBV & HCV are inactivated in 1 minute)

INFECTIOUS MATERIALS DISPOSAL: Any materials used in the cleaning process that may contain blood/body fluids are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

EFFICACY TESTS HAVE DEMONSTRATED THAT THIS PRODUCT IS AN EFFECTIVE BACTERICIDE & FUNGICIDE IN THE PRESENCE OF ORGANIC SOIL (5% BLOOD SERUM) FOR ULTRASOUND TRANSDUCERS & PROBES: Use only on external surfaces of ultrasound transducers and probes.

With a clean washcloth, clean gel from external surface of the ultrasound transducer or probe to be disinfected. Then thoroughly wet probe with wipe. Treated surface must remain wet for 5 minutes or refer to the contact times indicated in the chart. Use an additional wipe if needed to assure continuous 5 minute wet contact time. Let air dry.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

For surfaces that may come in contact with food, a potable water rinse is required.

For plastic and painted surfaces, spot test on an inconspicuous area before use. This product will not harm most surfaces: including acrylic, fiberglass and vinyl.

Not recommended for use on unpainted wood. Not for personal cleansing. This is NOT a baby wipe! Do not flush in toilet.

Fast, Effective & Convenient

Antibacterial Formula
Alcohols & Disinfectants

Alcohol & Bleach Free

Great On-the-Go

Cleaning

Lint-Free

Fresh Citrus

Scent



Disinfectant Wipes

All Purpose Cleaner

Cleans, Deodorizes, Disinfects & Sanitizes

Kills 99.9%

of Bacteria* in 15 seconds
of *Campylobacter jejuni* (ATCC 29428) in 1 minute
of *Listeria monocytogenes* (ATCC 35152) in 1 minute

Disinfects household, kitchen & bathroom bacteria!

Kills SARS-CoV-2 in

(the virus that causes COVID-19)



Active Ingredients

Alkyl (60% C14, 30% C16, 5% C12, 5% C18)
Dimethyl benzyl ammonium chloride..... 0.14% w/w
Alkyl (68% C12, 32% C14)
Dimethyl ethylbenzyl ammonium chloride... 0.14% w/w

DIN - 02514109

Lot #

Expiry Date:

NET CONTENTS

80 Wipes 19cm x 18.5cm
Single Wipe Size

* *Staphylococcus aureus* (Staph), *Klebsiella pneumoniae*, *Listeria monocytogenes*, *Escherichia coli* (E. coli) O157:H7
† *Pseudomonas aeruginosa*, *Staphylococcus aureus* (Staph), *Salmonella enterica* (Salmonella), *Streptococcus pyogenes* (Strep), *Escherichia coli* (E. coli) O157:H7, *Community Associated Methicillin Resistant Staphylococcus aureus* CA-MRSA (Genotype USA400), *Community Associated Methicillin Resistant Staphylococcus aureus* CA-MRSA (Genotype USA300), *Vancomycin resistant Enterococcus faecium* (VRE), *Methicillin resistant Staphylococcus aureus* (MRSA), *Vancomycin intermediate resistant Staphylococcus aureus* (VISA)

bullfrogproducts.ca

BullFrog Products
PO Box 77004
St. Albert, AB, T8N 6C1

Questions? Comments? 1-866-543-5276

CAUTION

See side panels for additional precautionary statements.

Read label before using.

KEEP OUT OF REACH OF CHILDREN

Disinfectant Wipes are non-porous, disposable wipes, pre-saturated with a quaternary disinfectant cleaner. **Disinfectant Wipes** is intended for use in:
• Hospitals & Healthcare Areas: Ambulances, Transport Vehicles, Anesthesia, CAT Labs, X-Ray, Surgery, Laboratory, Central Supply, Emergency Rooms, Operating Room, Recovery Room, CCU, ICU, Neonatal Intensive Care Units (NICU), Pediatrics, Pediatric Intensive Care Units (PICU), Home Health Care, Clinics, Dialysis Clinics, Hospices, Laundry Rooms, Newborn Nursery, Nursing Homes, Dental Office, Physician's Office, Doctor's Office, Optimal Office, Physiotherapy, Physical Therapy, Radiology, Rehabilitation
• Critical Care Areas: CCU, Emergency Rooms, ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU), Surgery, Ophthalmic Offices
• Veterinary Clinics, Animal Hospitals, Animal Life Science Laboratories, Kennels, Pet Animal Kennels, Breeding Establishments, Grooming Establishments, Pet Animal Quarters, Zoos, Pet shops, Tack shops
• Animal Holding Areas, Veterinary Waiting Rooms, Veterinary Examination Rooms, Veterinary X-Ray Rooms, Veterinary Office, Physiotherapy, Physical Therapy, Radiology, Rehabilitation
• Beauty Salons, Barber Shops, Correctional Facilities, Funeral Homes, Mortuaries
• Schools, Daycare Centers, Health Clubs, Gyms, Athletic Facilities, Airports, Transportation Terminals
• Grocery Stores, Shopping Centers, Delicatessens, Catering Facilities, Banquet Halls
• Household & Institutional Kitchens and/or Bathrooms

Disinfectant Wipes will clean & disinfect hard, non-porous, surfaces of: bed railings, cabinets, carts, chairs, counters, exam tables, gurneys, isotherms, IV poles, stethoscopes, stretchers, tables, telephones, toilet seats, toilets, bathtubs, bathroom fixtures, shower stalls, sinks, tiles, vanity tops, hospital beds, traction devices, MRI, CAT, scales, paddles, wheelchairs, ultrasound transducers & probes.

Disinfectant Wipes will clean & disinfect hard, non-porous surfaces of: grocery cart handles, grocery carts, shopping carts, grocery cart child seats.

Disinfectant Wipes will clean & disinfect hard, non-porous surfaces of: changing tables, cribs, diaper pails, trashcans, garbage cans.

Disinfectant Wipes will clean & disinfect hard, non-porous surfaces of: cabinets, counters, chairs, computers, desks, desktops, doorknobs, fax machines, garbage/trash cans, headsets, keyboards, railings, seats, sinks, tables, telephones, appliance exteriors, faucets, bathtubs, shower stalls, stove tops.

Disinfectant Wipes will clean & disinfect hard, non-porous surfaces of: gym equipment, exercise machines, wrestling mats, gymnastic equipment.

Disinfectant Wipes to clean & disinfect non-porous personal protective safety equipment, protective headgear, hard hats, half mask respirators, full face breathing apparatus, gas masks, goggles, spectacles, face shields, silicone rubber or PVC hearing protectors and vinyl covered ear muffs. Rinse all equipment that is prolonged contact with skin before reuse with clean warm water (about 120°F), and allow to air dry. (Caution: Cleaning at the recommended 120°F temperature will avoid overheating and distortion of the personal safety equipment that would necessitate replacement.)

Disinfectant Wipes is formulated to disinfect hard, non-porous, inanimate environmental surfaces such as floors, walls, metal surfaces, stainless steel surfaces, glazed porcelain, glazed ceramic tile, plastic surfaces, bathrooms, shower stalls, bathtubs, and cabinets.

VIRUCIDAL ACTIVITY		BACTERICIDAL ACTIVITY	
HIV-1 (associated with AIDS)	30 seconds	<i>Pseudomonas aeruginosa</i> (ATCC 15442)	2 minutes
Influenza A (H1N1) Virus (ATCC VR-1469)	1 minute	<i>Staphylococcus aureus</i> (ATCC 6538)	2 minutes
Influenza A2/05/57 (H2N2) (Charles River Laboratories A2/05/57)	2 minutes	<i>Salmonella enterica</i> (ATCC 10708)	2 minutes
Avian Influenza strain (H5N1)	30 seconds	<i>Escherichia coli</i> O157:H7 (ATCC 35150)	2 minutes
Hepatitis B Virus (HBV)	1 minute	<i>Streptococcus pyogenes</i>	5 minutes
Hepatitis C Virus (HCV)	1 minute	<i>Klebsiella pneumoniae</i> (New Delhi Duck Hepatitis B Virus (DHBV))	2 minutes
Duck Hepatitis B Virus (DHBV)	1 minute	<i>Metallo-Beta Lactamase</i> (NDM-1 Carbapenem Resistant)	2 minutes
Bovine Viral Diarrhea Virus (BVDV)	1 minute	ANTIBIOTIC-RESISTANT	
Human Coronavirus (ATCC VR-740)	2 minutes	<i>Methicillin resistant Staphylococcus aureus</i> (MRSA) (ATCC 35593)	2 minutes
SARS-Related Coronavirus 2 (SARS-CoV-2) (the virus that causes COVID-19) (NIA-WA1/2020-B-NR-02821)	30 seconds	<i>Vancomycin resistant Enterococcus faecalis</i> (VRE) (ATCC 51755)	5 minutes
Herpes Simplex Virus 1 (HSV-1)	2 minutes	<i>Vancomycin intermediate resistant Staphylococcus aureus</i> (VISA) (CDC HP 5836)	2 minutes
Herpes Simplex Virus 2 (HSV-2) (ATCC VR-734)	2 minutes	<i>Community Associated Methicillin resistant Staphylococcus aureus</i> (CA-MRSA) (NSR 129)	2 minutes
Rotavirus (ATCC VR-2018)	1 minute	Genotype USA400	2 minutes
Norovirus (Norwalk virus)	2 minutes	<i>Community Associated Methicillin resistant Staphylococcus aureus</i> (CA-MRSA) (NSR 394)	2 minutes
Feline Calicivirus (ATCC VR-782)	2 minutes	Genotype USA300	2 minutes
Phine Chinovirus Type 14 (ATCC VR-284)	2 minutes	FUNGICIDAL ACTIVITY	
Rhinovirus Type 39 (ATCC VR-340)	2 minutes	<i>Trichophyton interdigitale</i>	5 minutes
Non-Food Contact Hard Surface Sanitization			
<i>Escherichia coli</i> O157:H7 (ATCC 35150)	15 seconds		
<i>Staphylococcus aureus</i> (ATCC 6538)	15 seconds		
<i>Klebsiella pneumoniae</i> (ATCC 4352)	60 seconds		
<i>Campylobacter jejuni</i> (ATCC 29428)	60 seconds		
<i>Listeria monocytogenes</i> (ATCC 35152)	60 seconds		

STORAGE & DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

STORAGE: Store in an area inaccessible to children. Store at room temperature.

DISPOSAL: Dispose of wipe in trash after use. Do not flush. Tightly close lid between uses to retain moisture. For information on disposal of unused, unwanted product, contact the manufacturer or the appropriate municipal/provincial/territorial agency or product stewardship organization. Rinse the emptied container thoroughly prior to disposal. Dispose of the empty container in accordance with municipal/provincial/territorial requirements. Offer for recycling, if available.

PRECAUTIONARY STATEMENTS

Hazards to Humans & Domestic Animals.



CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling.

FIRST AID

Have the product container or label & Drug Identification Number (DIN) with you when calling a poison control center or doctor or going for treatment.

FOR EYE CONTACT: Hold eyes open & rinse slowly & gently with water for 15 - 20 minutes. Remove contact lenses if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.



French/Français Label



Manufactured by:
AxSys Direct MFG.
Edmonton, AB T6E 5Y8
axsysdirect.com